



Laboratory Medicine in the Era of Disruptive Technology

LMCE 2017 & KSLM 58th Annual Meeting

October 18-20, 2017
Grand Walkerhill Seoul, Korea
www.lmce-kslm.org

BD

Product Highlights

BD Phoenix M50



The BD Phoenix ID/AST solution incorporates the use of an oxidation reduction indicator, turbidometric growth detection, full on-panel antimicrobial concentrations with delayed resistance features and uses integrated software that can assist with antimicrobial stewardship. These ensure a complete solution for rapid, accurate and reliable AST results

BD FACSLyric™



Consistency you can rely on

The BD FACSLyric™ system combines automation, simplicity and reliability to help your laboratory be more productive. This powerful flow cytometer includes user friendly BD FACSuite™ Clinical software to make daily tasks easy and fast. In conjunction with a broad range of reagents and BD's history of service and quality, the BD FACSLyric system offers a truly integrated solution and provides clinicians with accurate, reliable and repeatable results – test to test, instrument to instrument and site to site, regardless of the assay complexity or operator experience level.