

LMCE 2017 & KSLM 58th Annual Meeting

October 18-20, 2017 Grand Walkerhill Seoul, Korea www.lmce-kslm.org

Thermo Fisher Scientific

Product Highlights

Ion Torrent Ion S5 Systems (Display)



Ion S5 Sequencer and Ion S5 XL Sequencer are DNA sequencing devices enable DNA analysis by measuring hydrogen ions released during nucleotide bindings via semi-conductor chips. This device is used to analyze DNA sequences at various specific regions among multiple human genes.

Phadia 250



Fully automated, Phadia 250 is ideal for the medium-sized laboratories delivering 400 - 2,000 results per week.

- An allergen panel detecting sensitization to 650 allergens and 90 allergen components
- A broad autoimmunity menu covering clinically relevant autoimmunity markers for evaluation of the most common autoimmunity diseases, e.g., rheumatoid arthritis, celiac disease and connective tissue diseases
- All tests can be measured from one sample
- User-configurable reflex testing functionality

