



KM Healthcare Corp.

Product Highlights

Temperature Monitoring Transport System



Self-sufficient and controlled refrigeration for long duration transport of refrigerated or frozen products. : The Coldway solution is to store dry ice or eutectic plates separately inside the container, in order to keep it away from the transported products. Thanks to its fan cooling system, inside temperature can be set between -30°C and $+8^{\circ}\text{C}$ with a precision band of $\pm 2^{\circ}\text{C}$. The Crypak's data logger generates a comprehensive PDF report without the need for proprietary software, making data retrieval faster and more convenient. (Temperature, Humidity)